

## Endoscopy Mask

Allows simultaneous ventilation and endoscopy of the patient to increase patient comfort and the working quality of the doctor.



- **Fiberoptic Intubation**

The ET Tube can be railroaded orally or nasally along the bronchoscope through the flexible membrane. After intubation the membrane is pushed inside which allows easy removal of the mask.

- **Bronchoscopy**

The membrane guarantees an ideal seal and therefore an efficient ventilation before, during and after bronchoscopy without time pressure.

- **Gastroenterology**

- **Trans Esophageal Echocardiography**

The use of Entonox (50% N<sub>2</sub>O, 50% O<sub>2</sub>) is ideal for minor procedures where pain is predictable and of short duration. Entonox can be self administered from the conscious patient via the mask. It has a rapid analgesic effect and therefore allows a pain free insertion of the Endoscope through the hole in the membrane. Entonox wears off quickly without post operative problems.

Product Code	Description	QTY
30-40-000	Endoscopy Mask for single use, Material: PVC Size 0, Newborn 2.0 mm	BOX 6
30-40-111	Endoscopy Mask for single use, Material: PVC Size 1, Infant 2.0 mm	BOX 6
30-40-333	Endoscopy Mask for single use, Material: PVC Size 3, Child 3.0 mm	BOX 6
30-40-335	Endoscopy Mask for single use, Material: PVC Size 3, Child 5.0 mm	BOX 6
30-40-555	Endoscopy Mask for single use, Material: PVC Size 5, Adult 5.0 mm	BOX 6
30-40-777	Endoscopy Mask for single use, Material: PVC Size 5, Adult 10.0 mm	BOX 6
30-40-400-1	Bronchoscope Airway for single use, sterile, Material: EVA Size 0, 36 mm	BOX 10
30-40-420-1	Bronchoscope Airway for single use, sterile, Material: EVA Size 2, 65 mm	BOX 10
30-40-440-1	Bronchoscope Airway for single use, sterile, Material: EVA Size 4, 85 mm	BOX 10
30-40-266	Hook ring, Ø 32 mm fits for Endoscopy Masks, for single use, Material: PA	BOX 10
35-30-255	Fixation tape for single use, Material: Isoprene, Adult	BOX 10
35-70-155	Fixation tape autoclavable, Material: Silicone, Child	BOX 10
35-70-255	Fixation tape autoclavable, Material: Silicone, Adult	BOX 10
35-40-622	Fixation tape for single use, Material: Velours / Cotton, For mask # 1	BOX 10
35-40-633	Fixation tape for single use, Material: Velours / Cotton, For mask # 3, 5	BOX 10