

Rust/Stain Removers & Descalers



Surgistain® Item No. 34572

Removes Rust, Stains and Spotting from Stainless Steel Surgical Instruments

Routine use will enhance the life and efficiency of stainless steel instruments

- Safe, efficient and quick revitalizing solution that removes rust, stains, spotting, hard water scale and mineral deposits from stainless steel surgical instruments
- Perfect for use on stainless steel trays, basins, and carts
- Routine use will enhance the life and efficiency of surgical instruments helping to loosen stiff joints and locks
- Excellent for use in ultrasonic machines
- Perfect for use in Operating Rooms, Outpatient Surgery, Ophthalmology, Dental and other departments for instrument restoration
- Free rinsing; non-aerosol; environmentally safe



F*A*C*T*S*® Item No. 34533

Foaming Autoclave Cleaning Treatment System

- Removes all rust, stains and corrosion from walls of autoclaves and from stainless steel carts, and trays
- Also removes scale, alkaline, and mineral deposits
- Foaming spray remains on walls longer for better penetration
- Eliminates all noxious odors
- Restores original finish to autoclaves, carts, and trays
- Keeps carts and trays looking new
- Non-corrosive; will not harm quality stainless steel



A*C*T*S*® Item No. 34501

Non-foaming formula of
F*A*C*T*S*®

Please visit www.ruhof.com for additional product information and current SDS.

Rust/Stain Removers & Descalers



Delyme® Item No. 34507

Lime and Mineral Scale Remover

Specifically designed for use in all automatic washers, washer sterilizers, cart washers and tunnel washers

- Removes hard water scale, mineral scale, alkaline mineral deposits, rust, stains and organic waste from the inside walls of automatic washers
- Regular use will prevent the build up of lime and mineral scale
- Helps maintain the condition of your washer disinfectant equipment
- Improves the cleaning process

Please visit www.ruhof.com for additional product information and current SDS.

Lubricants



Premixslip® Item No. 34552

Instrument Lubricant and Rust Inhibitor

Thoroughly steam penetrable and steam sterilizable

- High-grade mineral oil lubricant and rust inhibitor.
- Prevents rust, spots, stains and corrosion thus prolonging the life of the instrument.
- Prevents damage due to the misuse of stiff or grinding instruments thus reducing repair and replacement costs.
- Steam penetrable, steam sterilizable and ETO sterilizable.
- Non-oily, non-sticky, non-toxic and silicone free.



Surgislip® Item No. 34571

Concentrated Instrument Lubricant and Rust Inhibitor for Use in Washers



Premixslip® AW Clear Item No. 34555

Instrument Lubricant and Rust Inhibitor for Use In All Washers

Prohibits rusting, staining and spotting on all stainless steel instruments while adding a protective coating to prevent corrosion

- Perfect for use in all washers
- Prevents rust, spots, stains and corrosion thus prolonging the life of the instrument.
- Prevents damage due to misuse of stiff or grinding instruments reducing repair and replacement costs.
- Steam penetrable, steam sterilizable and ETO sterilizable.
- Non-oily, non-sticky and non-toxic; Free of silicone and mineral oil.

Please visit www.ruhof.com for additional product information and current SDS.

Cart & Tray Cleaners



ISO-GON® Item No. 34537

Versatile Cleaner, Descaler, Decontaminant, and Radioactive Contaminant Remover

Designed for use in washer-disinfectors, tunnel washers and cart washers

- One-step cleaning and scale prevention
- Mildly acidic (pH 5-5.5) yet does not require a neutralizer.
- Extremely fast wetting, while maintaining high detergency and cleaning performance in hard, soft, cold and hot water
- Bio-degradable
- Instantly and constantly active (sequestering)
- Free rinsing, leaves no residue or pyrogens



Rinse Aid® Item No. 34563

Prevents Spotting on Surfaces of all Surgical Instruments and Glassware

Designed for all automatic washers, washer sterilizers, cart washers and tunnel washers

- Designed for all automatic washers
- Prevents adhesion of mineral deposits
- Prevents spotting
- Assists in drying
- Leaves no residue

Please visit www.ruhof.com for additional product information and current SDS.