





HADDENHAM®

Comfiwave®




not just a night-time compression garment

What is *Comfiwave*[®]?

Comfiwave is a patented compression garment, designed to complement the treatment a patient receives during waking hours.

Comfiwave gives a therapeutic pressure of 15-21mmHg.

Comfiwave is available in below knee and thigh high styles for the lower limb; and with or without mitten for the upper limb.



More than just a night garment...

Easy donning and doffing gives Comfiwave an innovative versatility, providing both night-time compression and a comfortable solution when resting in the daytime.

Made from 82% Cotton;
10% Polyurethane; and 8% Polyester.

Features and benefits



Easy to don and doff

Patented, high-elastic knit makes Comfiwave easy to put on and take off



Made from a soft, undyed cotton

100% cotton against the skin and seamless design ensures maximum comfort



Breathable fabric

Comfiwave's wave patterning gives it high air permeability for coolness and stretch



Ambidextrous

Designed to be used on either the left or right limb

Indications



Provides safe, night-time compression



Can be used in palliative care



Can be worn under easywrap or bandaging




Compression to drift-off in

Comfiwave®










































































Lower limb

15-21 mmHg

Circumferences cm		S	M	L	XL
g	THIGH TOP	40-48	46-56	54-65	62-75
c	WIDEST CALF	25-32	30-38	36-46	44-56
b	ANKLE	17-22	20-26	22-30	27-35



Lengths cm		SHORT	REGULAR	LONG
A-d	BELOW KNEE	32-37	37-42	42-48
A-g	THIGH HIGH	64-70	70-77	77-84


		SHORT	REGULAR	LONG
 BELOW KNEE \$237.00	S	CW1-AD2S  	CW1-AD2R  	CW1-AD2L    
	M	CW1-AD3S  	CW1-AD3R   	CW1-AD3L    
	L	CW1-AD4S  	CW1-AD4R   	CW1-AD4L    
	XL	CW1-AD5S  	CW1-AD5R   	CW1-AD5L    
 THIGH HIGH \$363.00	S	CW1-AG2S  	CW1-AG2R   	CW1-AG2L    
	M	CW1-AG3S  	CW1-AG3R   	CW1-AG3L    
	L	CW1-AG4S  	CW1-AG4R   	CW1-AG4L    
	XL	CW1-AG5S  	CW1-AG5R   	CW1-AG5L    

Comfiwave®

































Upper limb

15-21 mmHg

Circumferences cm		S	M	L
g	AXILLA	25-30	30-40	40-50
d	MID-POINT OF FOREARM	23-26	26-30	30-38
c	WRIST	13-18	18-22	22-28
a	PALM AT BASE OF FINGERS	16-22	18-24	20-26



Lengths cm		SHORT	REGULAR
c-g	WRIST TO AXILLA	36-42	40-46

		SHORT	REGULAR
 <p>SLEEVE \$205.00</p>	S	CW1-CG2S  	CW1-CG2R   
	M	CW1-CG3S  	CW1-CG3R   
	L	CW1-CG4S  	CW1-CG4R   
 <p>SLEEVE + MITTEN \$337.00</p>	S	CW1-AH2S  	CW1-AH2R   
	M	CW1-AH3S  	CW1-AH3R   
	L	CW1-AH4S  	CW1-AH4R   

CWLF-001 ver.01

Tel: 0800 333103

Fax: 09 622 1234

Email: orders@jackson-allison.co.nz

Website: www.jackson-allison.co.nz

**JACKSON ALLISON MEDICAL & SURGICAL LTD.
56 Lunn Avenue, Mt Wellington, Auckland 1072
PO BOX 4261, Shortland Street, Auckland 1140
New Zealand**