

**JACKSON
ALLISON**

Aspen Surgical

Anti-Fatigue Mat



**Aspen
Surgical.**

Intended to act as an ergonomic product that reduces fatigue of nurses and surgeons experienced from long hours of standing in one location.

If a caregiver is standing in the same position for more than two hours continuously or more than 30% of the work day then they should use fatigue reducing techniques (i.e. alternate propping one foot on footstool, anti-fatigue mats, sit-stand stool, supportive footwear).

2012 AORN Perioperative Standards and Recommended Practices pg. 698.

FEATURES

- ✓ Intended to act as an ergonomic product that reduces fatigue of nurses and surgeons experienced from long hours of standing in one location
- ✓ Uses a dual foam and gel technology to provide comfort under foot
- ✓ Top surface provides necessary friction level with damp or dry surgical booties
- ✓ Bottom surface is made of a high-traction material to provide slip resistance
- ✓ Sides are beveled to prevent tripping and it is a bright red color for easy visibility
- ✓ Compatible with hospital grade quaternary cleaners
- ✓ Water Resistant
- ✓ Reusable for 6-12 months
- ✓ Non-Sterile
- ✓ Not made with natural rubber latex




Reducing fatigue of nurses and surgeons during procedures.

Anti-microbial and anti-bacterial
Cal 117 compliant

Gel core

High-performance foam

High-traction bottom surface
certified by the National
Floor Safety Institute



ORDERING INFORMATION

Product Number	Description	QTY
81005	Anti-Fatigue Mat Red 13" x 15"	BOX OF 5
82001	Anti-Fatigue Mat Red 18" x 24"	BOX OF 1
82005	Anti-Fatigue Mat Red 18" x 24"	BOX OF 5
83002	Anti-Fatigue Mat Red 20" x 48"	BOX OF 2
84003	Anti-Fatigue Mat Red 20" x 32"	BOX OF 3

**JACKSON
ALLISON**

Jackson Allison Medical & Surgical Ltd
56 Lunn Ave, Mt Wellington, Auckland 1072
PO Box 4261, Shortland St, Auckland 114

P: + 0800 333103 F: + 09 622 1234
Email: orders@jackson-allison.co.nz
Website: www.jackson-allison.co.nz