

Compression Stocking/Garment Applicator

**Works with both open
and closed toe stockings**

**Fits Compression
Garments to Upper
and Lower Limbs**

Significantly Places Less Stress and Pain on the Injury



ezy-as


Australian Made and Owned



EZY-AS™ "The Aussie Solution"

WHAT IS IT?...Ezy-As™ is a unique Australian Made product designed to assist users with the application of compression stockings, tubular bandages and arm garments.

HOW DOES IT WORK?...The compression stocking or tubular bandage is turned inside-out over the applicator thus distributing the elastic forces of the garment onto the smooth outer surface of the device. Once the garment is correctly positioned on the applicator the user then transfers the compression stocking, tubular bandage or arm garment onto the limb. For more information including: video demonstrations, instruction guide, testimonials... please visit us at: www.ezyasabc.com



Available in three colour coded sizes for easy identification




Robust, durable and easy to clean

Practical for both open and closed toe stockings

Compact, lightweight and travel friendly

ORDER CODE AND SIZE INFORMATION



Code*	Size	Colour	Measurement A or B
SM259	Small	Red	cm (Inches)
SM260	Medium	Yellow	28-35 (11-13.5)
SM261	Large	Blue	32-40 (12.5-15.75)
SM262	Handle	Black	38-47 (15-18.5)
			Fits all sizes



Optional Handle sold separately

* Codes used for orders placed with Sutherland Medical Pty Limited



EZY-AS ABC Pty Ltd
ACN 120 626 486 ABN 67 120 626 486

For enquiries and orders:



Sutherland Medical Pty Limited

8 Kingston Town Close, Oakleigh Victoria 3166
P.O. Box 1194 Huntingdale Victoria 3166
Phone: 1300 664 027 Fax: 1300 664 028
Email: orders@sutherlandmedical.com.au
www.sutherlandmedical.com.au